

## **Protocol for cDNA synthesis and qRT-PCR**

### **cDNA Synthesis**

Superscript III 1<sup>st</sup> Strand Synthesis Kit (Invitrogen, Cat #18080-051)

We typically use 1µg of Total RNA per 20µl cDNA reaction, but you may be able to use less.

Initial Primer/RNA Mix:

0.5 µl random hexamer primers

0.5 µl oligo(dT) primers

1µl dNTP mix

8ul of RNA (1µg) + DEPC H<sub>2</sub>O

Heat to 95°C, 5 min

Buffer/Enzyme mix (assemble in this order):

2µl 10x cDNA synthesis buffer

4µl MgCl<sub>2</sub>

2 µl 0.1 M DTT

1µ RNaseOut

1µl Superscript III

Mix Buffer/enzyme components and add to Primer/RNA mix in a 0.2 ml microfuge tube.

Place the tube in the PCR machine programmed as follows:

25°C, 10 min, 50°C, 50 min, 85°C, 5 min.

Add: 1 ul RNase H, mix, incubate at 37°C, 20 min, 95°C, 10 min (an additional step to the Invitrogen protocol to remove any components that may interfere with the PCR (Ambion website). Store cDNA at -20°C

### **Real Time-PCR**

We typically dilute the cDNA reaction 1:10 with DEPC water and use 5µl as template for a 20µl RT-PCR reaction.

1. Pipette 5µl of diluted cDNA into all wells.
2. Pipette 5µl of each primer pair mixture into triplicate wells.  
Note: if the experiment requires multiple plates, a control primer (housekeeping gene) is run in triplicate on each plate to control for plate to plate or machine to machine variations.
3. Pipette 10µl of Bio-Rad iQ™ SYBR Green supermix (Cat. #170-1882) into all wells for a total reaction volume of 20µl.
4. Seal plate with optical film. Centrifuge PCR plate for 2 minutes at ~2500rpm to get all liquids to the bottom of the wells.

**RT-PCR Amplification Protocol** (Bio-Rad myIQ™ or iCycler iQ™ machines):

1. Begin with a 3 min step at 95°C

2. 40 cycles with a step of 10 sec at 95°C, 45 sec at 65°C, and 20 sec at 78°C (Data acquisition step).
3. One step of 1 min at 95°C.
4. One step of 1 min at 55°C
5. 80 cycles of 10 sec each starting at 55°C with a 0.5°C increment at each step up to 95°C (Melting point step to check for primer dimer formation).
6. For iCycler machine: use the FAM-490 filter for SYBR green detection.

**Plate Setup** (Bio-Rad myIQ™ or iCycler iQ™ software):

1. Label triplicate wells starting with the first three wells being labeled 1, next three labeled 2, etc. up to 32 which is the maximum number of triplicate samples that can be tested on a single 96 well PCR plate.
2. Change standard quantity of all wells to 1.